# Luxatemp

## Automix Solar

## Inhalt | Contents

|                                    |    |
|------------------------------------|----|
| Gebrauchsinformation - DEUTSCH     | 4  |
| Instructions for use - ENGLISH     | 11 |
| Mode d'emploi - FRANÇAIS           | 18 |
| Istruzioni d'uso - ITALIANO        | 26 |
| Instrucciones de empleo - ESPAÑOL  | 33 |
| Instruções de uso - PORTUGUÊS      | 41 |
| Gebruiksaanwijzing - NEDERLANDS    | 48 |
| Brugsanvisning - DANSK             | 55 |
| Användarinstruktioner - SVENSKA    | 62 |
| Instrukcja użycia - POLSKI         | 69 |
| Информация по применению - РУССКИЙ | 76 |

## Produktbeschreibung

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### Product description

---

Luxatemp-Automix Solar is a light-curing composite for the fabrication of temporary crowns and bridges, inlays, onlays and veneers. The 2-component material is based on multifunctional methacrylates and is automatically mixed.

### Indications

---

- Fabrication of temporary crowns, bridges, inlays, onlays, partial crowns and veneers.
- Fabrication of long-term temporaries.

### Contraindications

---

Do not use the material in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

### Basic safety instructions

---

- Only for use by dental professionals!
- Keep away from children.

- Avoid contact with skin and eyes. In the event of inadvertent contact, rinse thoroughly and immediately with plenty of water and consult a physician if necessary.
- Wear a face mask and goggles when finishing the temporary, otherwise sanding dust may get into the eyes and respiratory tract.

## Application system

- Automix cartridge: see “Using the Automix cartridge” on page 14

**Note:** When using a cartridge for the first time, express a pea-sized amount of material and discard it.

## Timing

|               |                        |
|---------------|------------------------|
| 0:00–1:10 min | Insertion in the mouth |
| 2:00–4:00 min | Removal from the mouth |

**Note:** The times provided apply to a room temperature of 23 °C and to a normal relative air

humidity of 50%. Higher temperatures shorten, lower temperatures increase those times.

| Light curing       |               |
|--------------------|---------------|
| 20 sec per element | Hand lamp     |
| ≥ 1:00 min         | Light furnace |

## Recommended use

1. Before the crown or bridge preparation or extraction of teeth, take an impression using alginate, silicone or a thermoplastic impression material.

Alternatively, the temporary can be made using a laboratory-prepared thermoforming tray or, in the case of single crowns, with a strip crown. Strictly follow the appropriate manufacturer’s instructions!

**Note:** Using a silicone impression, undercuts have to be blocked out and, if necessary, pressure relief channels may have to be cut in. To create a connection with the abutment teeth in the

case of gaps in the molar region, cut a channel into the impression.

2. Dry the prepared teeth and block out any undercuts in the preparation. Then lightly lubricate the prepared teeth and surrounding tissue plus any composite core build-ups using Vaseline or a similar separating medium.
3. Place the Luxatemp-Automix Solar cartridge in the Automix-Dispenser (see “Using the Automix cartridge” on page 14).

**Note:** To avoid bubbles, it is important that the end of the mixing tip always remains immersed in the material during application.

4. Fill the mould or the thermoforming tray with Luxatemp-Automix Solar. First, apply the material to the occlusal surfaces of the mould or thermoforming tray and fill up to the gingival areas, only slightly overfilling in the process.

5. Place the impression or thermoforming tray, filled with Luxatemp-Automix Solar, over the prepared teeth within 1:10 minutes after start of mixing, using moderate pressure. Fix in place.

**Note:** The setting reaction has to be monitored intra-orally (e.g. using a scaler), because the curing is significantly affected by the temperature in the mouth and the temporary can only be removed during its elastic phase.

6. Remove the temporary from the mouth during the elastic phase, i.e. from 2:00 to 4:00 minutes after start of mixing.

**Note:** The temporary restoration must be light-cured in order to achieve definitive hardness.

7. Light curing: Expose the temporary restoration to a halogen lamp, 20 sec per element, or cure in a light furnace for at least 1:00 min.
8. Before finishing, remove the oxygen inhibition layer from the surface of the temporary (e.g. with alcohol).

**Caution!** When finishing the temporary, sanding dust may get into the eyes and respiratory tract.

- ▶ Wear face mask and goggles.
- ▶ Do not inhale the sanded dust and use an evacuator.

9. Refine the temporary using cross-toothed grinders, flexible discs, etc.
10. Check occlusion and grind if necessary using suitable instruments.
11. If necessary, the surface of the temporary may be glazed with a light-curing high-gloss varnish (e.g. Luxatemp-Glaze & Bond / LuxaGlaze in the USA) or polished manually with a rubber polisher or goat's hair wheel.
12. Clean and dry the prepared tooth and thoroughly remove the separating medium.
13. Cement the temporary.

**Note:** Any temporary cement can be used for cementing the temporary. DMG recommends using automatically mixing luting cements, e.g. TempoCemNE (non-eugenol zinc oxide cement).

Cements containing eugenol can inhibit the polymerization of resin-based luting cements during the permanent cementation.

14. Check occlusion and adjust accordingly.

## Customizing and repair of temporaries

---

**Note:** When customizing or repairing a freshly-made temporary it does not need to be roughened and a bonding agent is not required (proceed with step 2).

**Caution!** When working on the temporary, sanding dust may get into the eyes and respiratory tract.

- ▶ Wear face mask and goggles.
- ▶ Do not inhale the sanded dust and use an evacuator.


1. Repair of an older temporary: Roughen the temporary slightly using a grinder or a sand-blaster at the appropriate points. Then apply a bonding agent (e.g. Luxatemp-Glaze & Bond) according to manufacturer's instructions.

2. Apply Luxatemp-Automix Solar to all surfaces.
3. Hold the parts together for 3 minutes.
4. After 5 minutes, the appropriate points can be finished as usual.

**Note:** Alternatively, the temporary restoration may also be customised or repaired using the free-flowing and light-curing LuxaFlow Star composite or Luxatemp Star. Observe corresponding directions for use!

## Using the Automix cartridge

---

### Insert cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C], insert the cartridge and lock it by pushing down the plastic lever [C].

**Note:** Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

### Attach mixing tip

1. Remove the cartridge cap [D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

**Note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it into place.

## Apply material

- ▶ When using a cartridge for the first time, express a pea-sized amount of material and discard it.
- ▶ The material is mixed in the mixing tip by pulling the trigger [F] of the Automix-Dispenser and it is then applied directly.

**Note:** Leave the used mixing tip on the cartridge after use. It serves as a cap.

## Remove cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] completely.
2. Lift up the plastic lever [C] and remove the cartridge.

## Disinfection of the Automix-Dispenser

The Automix-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends

using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser can also be autoclaved.

## Side effects

There are no known systemic side effects to date.

## Composition

Glass filler in a matrix of multifunctional methacrylates; catalysts, stabilizers and additives. Free of methyl methacrylate and peroxides. Total filler volume: 44 w. % = 22 vol. % (0.02 to 2.5  $\mu\text{m}$ )

## Storage and shelf life

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F).
- ▶ Do not store in the refrigerator.
- ▶ Protect against exposure to direct sunlight.
- ▶ Leave the used mixing tip on the Automix cartridge after use. It serves as a cap.
- ▶ Do not use after the expiration date!

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH