LuxaCam PEEK Block

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English

Product description

LuxaCam PEEK block – tooth-colored, metal-free milling blocks for the fabrication of dental restorations. This material is not radiopaque.

Intended purpose

Fabrication of dental restorations using CAD/CAM milling machines:

- Fixed dental prostheses: Crown and bridge frameworks in the anterior and posterior region
- Removable dentures: telescopic primary and secondary crowns
- Individual implant abutments: Superstructures in the front and posterior region

Limitations of use

Bridge framework: with a maximum of two pontics between two abutments

Indications

- Tooth loss, edentulism
- Major tooth defects that cannot be restored by direct restoration, e.g. extensive carious lesions, traumatic injuries or extensive dental wear
- Endodontically treated teeth that are susceptible to fractures

Contraindications

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

People treated in the course of a dental procedure.

Intended users

Dentist, dental technician

Notes for use

Always process the material with the specified machine templates to avoid overheating and thus prevent damage to the material.

Recommended use

Processing the blocks

Restoration designs/design specifications

Design the dental structures for the restoration with suitable CAD software according to the manufacturer's information.

Creating the structure

Note: The blocks can be processed with all open CAD/CAM milling machines. Use a 5-axis milling machine to reach undercutting areas optimally.

Note: Milling/grinding processes and machine templates must be obtained from the machine manufacturers.

1. Secure the block into the CAD/CAM milling machine.
2. Start the milling process.

Removing the structure

Carefully remove the milled structure from the block with a separator or a milling tool.

Preparing the structure

1. Carry out the rough preparation, for example, with a carbide milling tool with FSQ toothing or a black ring cutter for titanium.

Caution! Material damage. When polishing, keep the pressure low to prevent smearing due to heat development.

2. Carry out the fine preparation and polishing with simple, crosscut HM milling tool.

Caution! Material damage. The following are unsuitable for cleaning: Steam, autoclaving, organic solvents and acids not expressly mentioned here, highly concentrated cleaning alcohols, chemical prosthesis cleaners, abrasive cleaning agents.

Veneering the structure

1. Individualize the structure as needed with an appropriate veneering resin, in accordance with the manufacturer's specifications.
2. Sandblast areas to be veneered (aluminum oxide, 50 – 110 μm , 2 – 3 bar).
3. Clean the restoration with alcohol.
4. Clean the structure with oil-free and water-free air.
5. For the veneering of PEEK, apply a suitable primer to the structure (e.g. Luxatemp-Glaze & Bond, DMG, in accordance with the manufacturer's specifications).

Luting the restoration

Temporary luting

Lute restorations with a suitable temporary cement or glass ionomer cement in accordance with manufacturer's instructions.

Permanent luting

Caution! Material damage. The following are unsuitable for cleaning: Steam, autoclaving, organic solvents and acids not expressly mentioned here, highly concentrated cleaning alcohols, chemical prosthesis cleaners, abrasive cleaning agents.

1. Clean the inside of the restoration with isopropyl alcohol.
2. Sandblast restoration (aluminum oxide, 50–110 μm , 2–3 bar).
3. Clean the restoration with oil-free and water-free air.
4. Apply Luxatemp-Glaze & Bond in accordance with the manufacturer's instructions and light-cure.
5. Lute the restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Interactions

Known reactions or interactions of the medical device with other materials already in place in the mouth should be taken into consideration by the dentist.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings / precautions

- For dental use only.
- Wear dust mask and goggles when milling, otherwise sanding dust may get into the eyes and respiratory tract.
- Keep out of the reach of children!
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse the eyes thoroughly with plenty of water and consult a physician if necessary.
- Keep out of the reach of children.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Performance characteristics of the device

Flexural strength ≥ 150 [MPa]
Flexural modulus ≥ 3.8 [GPa]

Storage / disposal

- Store in a dry place at room temperature (15–25 °C / 59–77 °F).
- Protect against exposure to direct sunlight!
- Do not use after the expiration date!
- Dispose of in accordance with national regulations.

Composition

Polyether ether ketone (approx. 80%), titanium dioxide (approx. 20%) with additives and pigments.

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